

Prothrombin time assay for measuring rivaroxaban plasma concentrations using calibrators and controls: results of a multicentre field trial

CPC047

Meyer Michel Samama^{1,2}, Genevieve Contant³, Theodore E Spiro⁴, Elisabeth Perzborn⁵, Léna Le Flem², Céline Guinet², Yves Gourmelin³, Jean Luc Martinoli³, for the Rivaroxaban Prothrombin Time Field Trial Laboratories

¹Hôtel-Dieu University Hospital, Paris, France; ²Biomnis Laboratories R&D, Ivry-sur-Seine, France; ³Diagnostica Stago SA, Gennevilliers, France; ⁴Bayer HealthCare Pharmaceuticals Inc., Montville, NJ, USA; ⁵Bayer HealthCare Pharmaceuticals, Wuppertal, Germany

Introduction

- Rivaroxaban – an oral, direct Factor Xa inhibitor¹ – is currently used in clinical practice for the prevention and treatment of thromboembolic disorders
- Routine coagulation monitoring is not required,² but a quantitative determination of rivaroxaban exposure might be useful in certain clinical circumstances (e.g. prior to urgent surgery)
- Because of its mode of action, rivaroxaban prolongs the prothrombin time (PT), but the results vary depending on the assay reagents; the international normalized ratio (INR) correction used for monitoring the vitamin K antagonists cannot be used for rivaroxaban^{3,4}

Objective

- To evaluate the suitability of the PT assay for the measurement of rivaroxaban plasma concentrations (ng/ml) using rivaroxaban calibrators and controls, and to assess the inter- and intra-laboratory precision of the measurements

Methods

- Participating laboratories in Europe and North America were provided with sets of rivaroxaban calibrators (0, 41, 219 and 430 ng/ml) and pooled human plasma controls containing 19, 160 and 643 ng/ml of rivaroxaban. The concentrations of rivaroxaban in the pooled human plasma controls were unknown to the participating laboratories
- Evaluations were carried out over 10 days by each laboratory using its own local PT reagent (Table 1) as well as a centrally provided PT reagent, STA[®] Neoplastine[®] CI Plus (Diagnostica Stago, Gennevilliers, France)
- Day-to-day precision and accuracy were evaluated by producing a calibration curve each day and by testing in duplicate the three pooled human plasma controls
- The control plasma samples were diluted with calibrator containing 0 ng/ml rivaroxaban when the values were greater than the calibration range

Results

- Local PT reagents:
 - A large inter-laboratory variation was seen when results were expressed in seconds; the coefficient of variation (CV) was 13.6–29.7%. Less variation was found when the results were expressed as rivaroxaban concentrations (ng/ml; CV 3.9–15.5%; undiluted samples), although over-estimation was observed (Figure 1; Table 2)
 - The intra-laboratory CV was 2.7–34.1% (for 19 ng/ml), 1.1–7.9% (160 ng/ml) and 1.1–9.6% (643 ng/ml)




Figure 1. Rivaroxaban concentrations in control plasma samples reported from individual laboratories using local prothrombin time reagents. Results are presented as median values \pm standard deviation (N=18).

Table 1. Local PT reagents and instruments used by the participating laboratories

Reagents	Central PT reagent	Local PT reagents					Total Local PT reagents	
		Diagnostica Stago			Siemens			Instrumentation Laboratory
		STA [®] Neoplastine [®] CI Plus	STA [®] Neoplastine [®] CI Plus	STA [®] Neoplastine [®] CI	Thromborel [®] S	Innovin [®]		
Diagnostica Stago	STA-R	7	3	1	1	1	6	
	STA Compact [®]	3	1			1	3	
Heinrich Amelung	KC 10	2	1			1	2	
Instrumentation Laboratory	ACL 1000	1		1			1	
	ACL TOP [®]	2				2	2	
Siemens	BCS	4	1		2	1	4	
Total		19	6	2	3	3	18	

PT, prothrombin time.

Central PT reagent:

- There was less inter-laboratory variation when the central PT reagent was used (CV 2.0–7.5%; undiluted samples; expressed as rivaroxaban concentrations) compared with local PT reagents; however, measured rivaroxaban concentrations were higher than the actual values (Figure 2; Table 2), as with local PT reagents
- The intra-laboratory CV was 4.5–19.3% (for 19 ng/ml), 1.2–8.3% (160 ng/ml), and 0.9–5.0% (643 ng/ml)
- The CV of the calibrators was 4.4–6.5% for the central PT reagent compared with 12.5–27.2% for local PT reagents (Table 2)




Figure 2. Rivaroxaban concentrations in control plasma samples reported from individual laboratories using the centrally provided prothrombin time reagent (STA[®] Neoplastine[®] CI Plus). Results are presented as median values \pm standard deviation (N=19).

Table 2. Inter-laboratory variations: comparison of results obtained with local PT reagents versus central PT reagent

	Actual rivaroxaban concentration	Local PT reagents		Central PT reagent (STA [®] Neoplastine [®] CI Plus)	
		Mean \pm SD (N=18)	CV, %	Mean \pm SD (N=19)	CV, %
Control samples, time (seconds)	19 ng/ml	12.6 \pm 1.7	13.6	13.8 \pm 0.8	5.9
	160 ng/ml	20.2 \pm 4.0	19.9	23.7 \pm 1.0	4.2
	643 ng/ml	44.7 \pm 13.3	29.7	57.2 \pm 2.5	4.4
	643 ng/ml*	23.9 \pm 5.9	24.9	28.5 \pm 1.8	6.4
Control samples, rivaroxaban concentration (ng/ml)	19 ng/ml	31 \pm 4	11.9	29 \pm 2	7.5
	160 ng/ml	197 \pm 31	15.5	186 \pm 4	2.2
	643 ng/ml	715 \pm 28	3.9	712 \pm 14	2.0
	643 ng/ml*	263 \pm 13	5.1	261 \pm 15	5.9
Calibrators, time (seconds)	0 ng/ml	11.0 \pm 1.38	12.5	11.8 \pm 0.77	6.5
	41 ng/ml	13.0 \pm 1.80	13.9	14.3 \pm 0.84	5.8
	219 ng/ml	22.1 \pm 5.09	23.0	26.5 \pm 1.20	4.6
	430 ng/ml	31.3 \pm 8.52	27.2	38.8 \pm 1.70	4.4

*Control samples were diluted twofold with calibrator (containing 0 ng/ml rivaroxaban). CV, coefficient of variation; PT, prothrombin time; SD, standard deviation.

Conclusions

- The results of this field trial suggest that it is feasible to measure rivaroxaban plasma concentrations (expressed in ng/ml) using the PT combined with rivaroxaban calibrators and controls, in contrast to the conventional INR, which cannot be used
- Owing to the variability of the measurements observed (in particular at low rivaroxaban plasma concentrations), more specific and sensitive methods (i.e. anti-Factor Xa chromogenic assays [please see poster CPC049]) are a better alternative when more precise measurements of rivaroxaban exposure are required
- Further validation of this method is required in clinical settings

References

- Perzborn E et al. *J Thromb Haemost* 2005;3:514–521.
- Bayer Pharma AG. Xarelto[®] (rivaroxaban) Summary of Product Characteristics. 2011. Available at: http://www.xarelto.com/html/downloads/Xarelto_Summary_of_Product_Characteristics_Dec2011.pdf. Accessed 2 March 2012.
- Samama MM et al. *Thromb Haemost* 2010;103:815–825.
- Lindhoff-Last E et al. *Thromb Haemost* 2010;32:673–679.

Disclosure of conflict of interest

This study was supported by Bayer HealthCare Pharmaceuticals and Janssen Research & Development, LLC (formerly Johnson & Johnson Pharmaceutical Research & Development, L.L.C.). Poster development was supported by Bayer HealthCare Pharmaceuticals and Janssen Research & Development, LLC. The data contained within this poster do not support or recommend the use of rivaroxaban in indications or countries in which it is not licensed.

Participants

DM Adcock, Englewood, CO, USA; JP Cambus, Toulouse, France; I Elalamy, Paris, France; G Gerotziakas, Paris, France; I Guoin-Thibault, Ivry-sur-Seine, France; J Harenberg, Mannheim, Germany; MH Horellou, Paris, France; M Jacquemin, Leuven, Belgium; K Jochmans, Brussels, Belgium; C le Courvoisier-Flaujac, Paris, France; C Legnani, Bologna, Italy; E Lindhoff-Last, Frankfurt, Germany; JL Martinoli, Gennevilliers, France; TL Ortel, Durham, NC, USA; G Palareti, Bologna, Italy; D Peetz, Mainz, Germany; MM Samama, Ivry-sur-Seine, France; A Schade, Cleveland, OH, USA; S Testa, Cremona, Italy; A Tripodi, Milan, Italy; AMHP van den Besselaar, Leiden, The Netherlands; J Weitz, Hamilton, Ontario, Canada; W Wijns, Brussels, Belgium.